

**BUY 3  
GET 1  
FREE!\***

## A World of White for Day and Night Specials to help you save

### Contrast**AM**<sup>®</sup>

#### In-office Whitening



ContrastAM is a 22% hydrogen peroxide gel designed for a one hour, in-office teeth whitening treatment. With a non-acidic pH, ContrastAM gel provides maximum efficacy and minimizes sensitivity. The dual barrel, 4:1 syringe delivery system also increases product stability. That means no refrigeration and a more convenient option for your practice.

ContrastAM is available in a Combo Kit that includes ContrastPM+ for at-home maintenance or a Power Pack, a great value that includes enough materials for 6 patient procedures.

### Contrast**PM+**<sup>®</sup>

#### At-home Whitening



ContrastPM+ at-home whitening gels are available in a variety of concentrations and kits to provide customized treatments for your patients. All ContrastPM+ formulas contain Potassium Nitrate, a desensitizing agent that helps deliver a more comfortable whitening experience. As a result, your patients are more likely to comply with treatment protocols and get optimal results.

ContrastPM+ 10%, 15% and 20% carbamide peroxide formulas are recommended for overnight wear (or 4-6 hours) in a custom tray. ContrastPM+ 36% has a recommended wear time of 15 minutes, twice-a-day and is ideal for patients who prefer not to sleep in bleaching trays or brux.

## ContrastAM In-Office Whitening

Item Description	Concentration	Applications
<b>Combo Kit with ContrastPM+ 15%</b> <i>Contents:</i> 1 Syringe ContrastAM whitening gel with tip (4.2ml) 1 Syringe FastDam® barrier material with tip 1 Syringe Contrast 1.1% fluoride* 1 Applicator brush 1 ContrastPM+ 15% touch-up kit 2 Sheets .040" EVA 1 Tray Case	22% HP	1 in-office procedure w/ touch-up
<b>Combo Kit with ContrastPM+ 36%</b> <i>Contents:</i> 1 Syringe ContrastAM whitening gel with tip (4.2ml) 1 Syringe FastDam® barrier material with tip 1 Syringe Contrast 1.1% fluoride* 1 Applicator brush 1 ContrastPM+ 36% touch-up kit 2 Sheets .040" EVA 1 Tray Case	22% HP	1 in-office procedure w/ touch-up
<b>Power Pack</b> <i>Contents:</i> 6 Syringes ContrastAM whitening gel with tips (4.2ml) 6 Syringes FastDam® barrier material with tips 6 Syringes Contrast 1.1% fluoride* 6 Applicator brushes	22% HP	6 in-office procedures



## ContrastPM+ At-Home Whitening

Item Description	Concentration	Applications
<b>Touch-Up Kit</b> <i>Contents:</i> 2 Filled Syringes (1.2ml ea.)	10% CP 15% CP 20% CP 36% CP <b>NEW!</b>	4 dual-arch
<b>Mini Kit</b> <i>Contents:</i> 4 Filled Syringes (1.2ml ea.)	10% CP 15% CP 20% CP 36% CP <b>NEW!</b>	8 dual-arch
<b>Full Patient Kit</b> <i>Contents:</i> 8 Filled Syringes (1.2ml ea.) 2 Sheets EVA 1 Tray Case	10% CP 15% CP 20% CP 36% CP <b>NEW!</b>	16 dual-arch



## Ancillary Products

Item Description	Contents
<b>FastDam® Soft-Tissue Isolation</b>	4 Filled Syringes (1ml ea.)
<b>Light-Curable Resin</b>	4 Filled Syringes (1.08ml ea.)
<b>EVA Tray Material</b>	12 Sheets (0.04in ea.)

